
L-715

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
INSTRUCTION MANUAL

LED COAXIAL EPISCOPIC LIGHT

Thank you for purchasing the HOZAN L-715 LED COAXIAL EPISCOPIC LIGHT. With proper care and handling, this fine instrument will provide years of trouble-free operation. Please read this entire instruction manual carefully before attempting to place this instrument in service. Please keep this instruction manual available for reference.

Identification of parts and contents

Specifications

LED	64 White LEDs
LED life	30000 Hrs * Varies to the used environment.
Color temp.	6500 K
Rating	12 Vdc (with Ac adapter)
Power consumption	3.8 W
Mounting support dia.	20 to 32 mm ϕ
Capacity	55.5 to 146.5 mm (with 20 mm ϕ post)
* Dimension from center of lens of inspection system to be applied to post	58 to 149 mm (with 30 mm ϕ post)
Dimension/window	32 x 32 mm

Dimensions	41(W) x 41(H) x 69(D) mm (Illumination unit) 60(W) x 34(H) x 108(D) mm (controller) 50(W) x 10(H) x 223(D) mm (bracket) *Without projections
Weight	135 g (Illumination unit) 130 g (Bracket)

Light brightness (at full power)

Distance from light source mm	Brightness lx
30	10300
60	4800
100	2300

Warning and caution symbols

These symbols are used throughout the instruction manual to alert the user to potential safety hazards as follows :

Warning ... Notice when incorrect handling could cause the user's death or serious injury.

Caution ... Notice when incorrect handling could cause injury to the user or material damage.

Even if the instructions do not have **Caution** mark, there are some possibilities for a serious situation. Follow the instructions.

Precautions

This device is a LED light for inspecting work. Do not use for any other purpose.


Caution

1. Use the provided ac adapter for power supply. Do not use other than this.
2. Do not apply a vibration or an impact to the device. Do not disassemble or modify the device.
3. The window of the top of the illumination unit is covered with glass. Avoid knocking it with a lens or such of inspection devices when work. The window at the bottom is not covered. Be careful not to hit the inner half-silvered mirror with specimens and such.
4. Do not touch the window glass and the inner half-silvered mirror with fingers directly. Gently wipe with soft cloth if fingerprints were left.
Do not look the ray directly. This could damage your eyes.


Operation

These products are used for this explanation.
L-835, L-815, L-521 and L-509

1 Install the bracket onto the post. Secure the height fix knob temporarily.

2 Put the illumination unit onto the bracket, and connect the controller and the ac adapter. The illumination unit must be placed turning its glassed face up as shown in the photo.
If the distance from the center of the lens to the post is less than 80 mm, place the illumination unit reversing the front and the rear (the camera screw faces the post).

* The illumination unit and the ac adapter can be connected directly without the controller, but LEDs could become more hot and shorten its life.

Ac adapter

Operation

- 3 Install the inspection devices, a camera and a lens and such. Adjust each position and fix to get focus.

Hint

Regarding the bracket position

The illumination unit must be parallel with the stage. Do not install the bracket tilting.

Regarding the illumination unit position (vertical)

Fix the illumination unit as near as possible to the lens. But, if the working distance (from the top of lens to the object) is long, it may need more.

As a guide; when fine images can not be gotten, e.g., an image is whitened and indistinct, with a working distance of exceeding 100 mm, adjust the position checking the images.

Regarding the illumination unit position (horizontal)

Fix the illumination unit aligning its center of the window and the center of the lens. Otherwise, fine image can not be gotten, e.g., it is out of focus.



Regarding the working distance

The working distance lengthen several mm by the characteristic of this device. Please focus after installing this device.

- 4 Use adjusting the brightness using the brightness adjust knob.


- 5 Switch off the power with the brightness adjust knob after use. Store the device disconnecting each connection.

The power is switched off at the dot position.

Move the fix screw corresponding the position of the illumination unit.

Illumination unit fix screw

Brightness adjust knob

Troubleshooting

On an occasion that inspection can not be done


If inspection can not be done enough using this device, please make sure of these points below.

Is there a trouble even without this device?	If so, make sure the conditions of the other devices than L-715.
Is not the illumination too bright or too dark?	Adjust the brightness operating the controller. (See Operation 4 on page 3.)
Is the illumination unit parallel with the stage?	Make sure the location of the bracket. (See Operation 3 on page 3.)
Is the illumination unit placed at the proper position?	Make sure the location of the illumination unit. (See Operation 3 on page 3.)

Mechanism of coaxial episcopic light

This device consists of mainly LEDs and a half-silvered mirror. The rays generated from the LEDs set at the side are reflected by the half-silvered mirror shine on the object which exist coaxially right under the camera sensor.

While the rays thrown on a horizontal face are reflected vertically and reach the sensor, the rays thrown on a rough face are reflected obliquely and do not reach the sensor. As a result, fine pictures with precise shade can be gotten.

HOZAN TOOL INDUSTRIAL CO.,LTD.

1-2-12 Saiwaicho, Naniwa-ku, Osaka 556-0021, Japan
Tel : +81-6-6567-3111 Fax : +81-6-6562-0024